



PH TEST PROCEDURE:

1. Fill test vial ½ full with solution to be tested.
2. Add 3 - 5 drops of pH test indicator to vial.
3. Tap the side of the test vial to mix contents.
4. Observe the color in the test vial.

Add pH DOWN (Acid) to lower pH.

Add pH UP (Base) to raise pH.

*The best pH is between 5.5 and 6.5.

pH **Test Indicator**

Wide Spectrum
pH 4.0 to 8.5
Super-Concentrated
Just 3-5 drops per test



**GENERAL HYDROPONICS
pH DOWN & pH UP**

are available in:

- liquid:
8oz, Quart, & 1 Gallon
- dry:
1 kg, 4lb, & 16lb

4.0
& lower

5.0

6.0

6.5

7.0

8.0

8.5
& higher

KEEP OUT OF REACH OF CHILDREN

STORE IN A COOL DARK PLACE

 **GH**
genhydro.com

Manufactured by GH Inc,
P.O. Box 1576 Sebastopol, CA 95473, USA



237 ml (8 fl oz)